




FLEXIBLE ENDOSCOPY ACCESSORIES



THE RZ FLEX LINE

TABLE OF CONTENTS

Biopsy Forceps.....	3-12
Foreign Body Forceps.....	13-18
Polyp Retriever, Scissors.....	19-26

Attention

Colour Code

- Ø 1.0 Schwarz / Black
- Ø 1.6 Weiss / White
- Ø 1.8 Blau / Blue
- Ø 2.1 Grün / Green
- Ø 2.3 Gelb / Yellow
- Ø 2.5 Schwarz / Black
- Ø 3.0 Rot / Red
- Ø 3.4 Orange / Orange



Long Oval Cup, fenestrated

Diameter	Length	Order-No.
1,0 mm	120 cm	310-901-120
1,0 mm	160 cm	310-901-160
1,0 mm	180 cm	310-901-180
1,0 mm	230 cm	310-901-230
1,4 mm	120 cm	310-902-120
1,4 mm	160 cm	310-902-160
1,4 mm	180 cm	310-902-180
1,4 mm	230 cm	310-902-230
1,6 mm	120 cm	310-903-120
1,6 mm	160 cm	310-903-160
1,6 mm	180 cm	310-903-180
1,6 mm	230 cm	310-903-230
1,8 mm	120 cm	310-904-120
1,8 mm	160 cm	310-904-160
1,8 mm	180 cm	310-904-180
1,8 mm	230 cm	310-904-230
2,1 mm	120 cm	310-905-120
2,1 mm	160 cm	310-905-160
2,1 mm	180 cm	310-905-180
2,1mm	230 cm	310-905-230
2,3 mm	120 cm	310-906-120
2,3 mm	160 cm	310-906-160
2,3 mm	180 cm	310-906-180
2,3 mm	230 cm	310-906-230

Standard delivery:

Metal tube distal tapered, metal handle.



On desire:

- ⇒ Metal tube Grounded, e.g. 310-000-230 S
- ⇒ Metal tube Isolated, e.g. 310-000-230 SI **NOTICE:** not for diameter 1.0 – 1.6 mm !
- ⇒ Plastic Handle, e.g. 310-000-230 P

ALL INSTRUMENTS ARE AUTOCLAVABLE
SPECIAL LENGTH ON REQUEST



Long Oval Cup, fenestrated, with needle

Diameter	Length	Order-No.
1,8 mm	120 cm	310-907-120
1,8 mm	160 cm	310-907-160
1,8 mm	180 cm	310-907-180
1,8 mm	230 cm	310-907-230
2,1 mm	120 cm	310-908-120
2,1 mm	160 cm	310-908-160
2,1 mm	180 cm	310-908-180
2,1 mm	230 cm	310-908-230
2,3 mm	120 cm	310-909-120
2,3 mm	160 cm	310-909-160
2,3 mm	180 cm	310-909-180
2,3 mm	230 cm	310-909-230

Standard delivery:

Metal tube distal tapered, metal handle.



On desire:

- ⇒ Metal tube Grounded, e.g. 310-000-230 S
- ⇒ Metal tube Isolated, e.g. 310-000-230 SI **NOTICE:** not for diameter 1.0 – 1.6 mm !
- ⇒ Plastic Handle, e.g. 310-000-230 P

ALL INSTRUMENTS ARE AUTOCLAVABLE
SPECIAL LENGTH ON REQUEST



Short Oval Cup, fenestrated

Diameter	Length	Order-No.
1,8 mm	120 cm	310-910-120
1,8 mm	160 cm	310-910-160
1,8 mm	180 cm	310-910-180
1,8 mm	230 cm	310-910-230
2,1 mm	120 cm	310-911-120
2,1 mm	160 cm	310-911-160
2,1 mm	180 cm	310-911-180
2,1 mm	230 cm	310-911-230
2,3 mm	120 cm	310-912-120
2,3 mm	160 cm	310-912-160
2,3 mm	180 cm	310-912-180
2,3 mm	230 cm	310-912-230
2,5 mm	120 cm	310-913-120
2,5 mm	160 cm	310-913-160
2,5 mm	180 cm	310-913-180
2,5 mm	230 cm	310-913-230
3,0 mm	180 cm	310-914-180
3,0 mm	230 cm	310-914-230
3,4 mm	180 cm	310-915-180
3,4 mm	230 cm	310-915-230

Standard delivery:

Metal tube distal tapered, metal handle.



On desire:

- ⇒ Metal tube Grounded, e.g. 310-000-230 S
- ⇒ Metal tube Isolated, e.g. 310-000-230 SI **NOTICE:** not for diameter 1.0 – 1.6 mm !
- ⇒ Plastic Handle, e.g. 310-000-230 P

ALL INSTRUMENTS ARE AUTOCLAVABLE
SPECIAL LENGTH ON REQUEST



Short Oval Cup, fenestrated, with needle

Diameter	Length	Order-No.
1,8 mm	120 cm	310-916-120
1,8 mm	160 cm	310-916-160
1,8 mm	180 cm	310-916-180
1,8 mm	230 cm	310-916-230
2,1 mm	120 cm	310-917-120
2,1 mm	160 cm	310-917-160
2,1 mm	180 cm	310-917-180
2,1 mm	230 cm	310-917-230
2,3 mm	120 cm	310-918-120
2,3 mm	160 cm	310-918-160
2,3 mm	180 cm	310-918-180
2,3 mm	230 cm	310-918-230
2,5 mm	120 cm	310-919-120
2,5 mm	160 cm	310-919-160
2,5 mm	180 cm	310-919-180
2,5 mm	230 cm	310-919-230
3,0 mm	180 cm	310-920-180
3,0 mm	230 cm	310-920-230
3,4 mm	180 cm	310-921-180
3,4 mm	230 cm	310-921-230

Standard delivery:

Metal tube distal tapered, metal handle.



On desire:

- ⇒ Metal tube Grounded, e.g. 310-000-230 S
- ⇒ Metal tube Isolated, e.g. 310-000-230 SI **NOTICE:** not for diameter 1.0 – 1.6 mm !
- ⇒ Plastic Handle, e.g. 310-000-230 P

ALL INSTRUMENTS ARE AUTOCLAVABLE
SPECIAL LENGTH ON REQUEST



Round Cup, fenestrated

Diameter	Length	Order-No.
2,3 mm	120 cm	310-930-120
2,3 mm	160 cm	310-930-160
2,3 mm	180 cm	310-930-180
2,3 mm	230 cm	310-930-230
3,4 mm	180 cm	310-931-180
3,4 mm	230 cm	310-931-230



Round Cup, fenestrated, with needle

Diameter	Length	Order-No.
2,3 mm	120 cm	310-932-120
2,3 mm	160 cm	310-932-160
2,3 mm	180 cm	310-932-180
2,3 mm	230 cm	310-932-230

Standard delivery:

Metal tube distal tapered, metal handle.



On desire:

- ⇒ Metal tube Grounded, e.g. 310-000-230 S
- ⇒ Metal tube Isolated, e.g. 310-000-230 SI **NOTICE:** not for diameter 1.0 – 1.6 mm !
- ⇒ Plastic Handle, e.g. 310-000-230 P

ALL INSTRUMENTS ARE AUTOCLAVABLE
SPECIAL LENGTH ON REQUEST



Round Cup, fenestrated 3.0 / 2.3 mm

Diameter	Diameter Cups	Length	Order-No.
2,3 mm	3,0mm	180 cm	310-933-180
2,3 mm	3,0mm	230 cm	310-933-230



Round Cup, fenestrated 3.0 / 2.3 mm, with needle

Diameter	Diameter Cups	Length	Order-No.
2,3 mm	3,0mm	180 cm	310-934-180
2,3 mm	3,0mm	230 cm	310-934-230

Standard delivery:

Metal tube distal tapered, metal handle.



On desire:

- ⇒ Metal tube Grounded, e.g. 310-000-230 S
- ⇒ Metal tube Isolated, e.g. 310-000-230 SI **NOTICE:** not for diameter 1.0 – 1.6 mm !
- ⇒ Plastic Handle, e.g. 310-000-230 P

ALL INSTRUMENTS ARE AUTOCLAVABLE
SPECIAL LENGTH ON REQUEST



Long Oval Cup, serrated

Diameter	Length	Order-No.
1,8 mm	120 cm	310-935-120
1,8 mm	160 cm	310-935-160
1,8 mm	180 cm	310-935-180
1,8 mm	230 cm	310-935-230
2,1 mm	120 cm	310-936-120
2,1 mm	160 cm	310-936-160
2,1 mm	180 cm	310-936-180
2,1 mm	230 cm	310-936-230
2,3 mm	120 cm	310-937-120
2,3 mm	160 cm	310-937-160
2,3 mm	180 cm	310-937-180
2,3 mm	230 cm	310-937-230



Long Oval Cup, serrated, with needle

Diameter	Length	Order-No.
1,8 mm	120 cm	310-938-120
1,8 mm	160 cm	310-938-160
1,8 mm	180 cm	310-938-180
1,8 mm	230 cm	310-938-230
2,3 mm	120 cm	310-939-120
2,3 mm	160 cm	310-939-160
2,3 mm	180 cm	310-939-180
2,3 mm	230 cm	310-939-230

Standard delivery:

Metal tube distal tapered, metal handle.



On desire:

- ⇒ Metal tube Grounded, e.g. 310-000-230 S
- ⇒ Metal tube Isolated, e.g. 310-000-230 SI **NOTICE:** not for diameter 1.0 – 1.6 mm !
- ⇒ Plastic Handle, e.g. 310-000-230 P

ALL INSTRUMENTS ARE AUTOCLAVABLE
SPECIAL LENGTH ON REQUEST



Serrated Jaw, v-shaped

Diameter	Length	Order-No.
1,8 mm	120 cm	310-940-120
1,8 mm	160 cm	310-940-160
1,8 mm	180 cm	310-940-180
1,8 mm	230 cm	310-940-230
2,3 mm	120 cm	310-941-120
2,3 mm	160 cm	310-941-160
2,3 mm	180 cm	310-941-180
2,3 mm	230 cm	310-941-230



Tip Jaw, v-shaped, 90° opened

Diameter	Length	Order-No.
1,8 mm	120 cm	310-942-120
1,8 mm	160 cm	310-942-160
1,8 mm	180 cm	310-942-180
1,8 mm	230 cm	310-942-230
2,3 mm	120 cm	310-943-120
2,3 mm	160 cm	310-943-160
2,3 mm	180 cm	310-943-180
2,3 mm	230 cm	310-943-230

Standard delivery:

Metal tube distal tapered, metal handle.



On desire:

- ⇒ Metal tube Grounded, e.g. 310-000-230 S
- ⇒ Metal tube Isolated, e.g. 310-000-230 SI **NOTICE:** not for diameter 1.0 – 1.6 mm !
- ⇒ Plastic Handle, e.g. 310-000-230 P

ALL INSTRUMENTS ARE AUTOCLAVABLE
SPECIAL LENGTH ON REQUEST



Biopsy Forceps, oval cups with fine teeth

Diameter	Length	Order-No.
1,8 mm	120 cm	310-944-120
1,8 mm	160 cm	310-944-160
1,8 mm	180 cm	310-944-180
1,8 mm	230 cm	310-944-230
2,3 mm	120 cm	310-945-120
2,3 mm	160 cm	310-945-160
2,3 mm	180 cm	310-945-180
2,3 mm	230 cm	310-945-230



Biopsy Forceps, oval cups with fine teeth and needle

Diameter	Length	Order-No.
1,8 mm	120 cm	310-946-120
1,8 mm	160 cm	310-946-160
1,8 mm	180 cm	310-946-180
1,8 mm	230 cm	310-946-230
2,3 mm	120 cm	310-947-120
2,3 mm	160 cm	310-947-160
2,3 mm	180 cm	310-947-180
2,3 mm	230 cm	310-947-230

Standard delivery:

Metal tube distal tapered, metal handle.



On desire:

- ⇒ Metal tube Grounded, e.g. 310-000-230 S
- ⇒ Metal tube Isolated, e.g. 310-000-230 SI **NOTICE:** not for diameter 1.0 – 1.6 mm !
- ⇒ Plastic Handle, e.g. 310-000-230 P

ALL INSTRUMENTS ARE AUTOCLAVABLE
SPECIAL LENGTH ON REQUEST



Biliary- and Pancreatic, Duct Forceps single action

Diameter	Length	Order-No.
2,3 mm	120 cm	310-948-120
2,3 mm	160 cm	310-948-160
2,3 mm	180 cm	310-948-180
2,3 mm	230 cm	310-948-230

Standard delivery:

Metal tube distal tapered, metal handle.



On desire:

- ⇒ Metal tube Grounded, e.g. 310-000-230 S
- ⇒ Metal tube Isolated, e.g. 310-000-230 SI **NOTICE:** not for diameter 1.0 – 1.6 mm !
- ⇒ Plastic Handle, e.g. 310-000-230 P

ALL INSTRUMENTS ARE AUTOCLAVABLE
SPECIAL LENGTH ON REQUEST



Alligator Forceps, cross serrated

Diameter	Length	Order-No.
2,3 mm	120 cm	310-950-120
2,3 mm	160 cm	310-950-160
2,3 mm	180 cm	310-950-180
2,3 mm	230 cm	310-950-230



Oval Cup, serrated

Diameter	Length	Order-No.
1,8 mm	220 cm	310-951-220
1,8 mm	180 cm	310-951-180

Standard delivery:

Metal tube distal tapered, metal handle.



On desire:

- ⇒ Metal tube Grounded, e.g. 310-000-230 S
- ⇒ Metal tube Isolated, e.g. 310-000-230 SI **NOTICE:** not for diameter 1.0 – 1.6 mm !
- ⇒ Plastic Handle, e.g. 310-000-230 P

ALL INSTRUMENTS ARE AUTOCLAVABLE
SPECIAL LENGTH ON REQUEST



Rat Tooth, 2:1 teeth

Diameter	Length	Order-No.
1,8 mm	120 cm	310-952-120
1,8 mm	160 cm	310-952-160
1,8 mm	180 cm	310-952-180
1,8 mm	230 cm	310-952-230
2,3 mm	120 cm	310-953-120
2,3 mm	160 cm	310-953-160
2,3 mm	180 cm	310-953-180
2,3 mm	230 cm	310-953-230
3,0mm	180cm	310-954-180
3,0mm	230cm	310-954-230

Standard delivery:

Metal tube distal tapered, metal handle.



On desire:

- ⇒ Metal tube Grounded, e.g. 310-000-230 S
- ⇒ Metal tube Isolated, e.g. 310-000-230 SI **NOTICE:** not for diameter 1.0 – 1.6 mm !
- ⇒ Plastic Handle, e.g. 310-000-230 P

ALL INSTRUMENTS ARE AUTOCLAVABLE
SPECIAL LENGTH ON REQUEST



Alligator Forceps

Diameter	Length	Order-No.
1,0 mm	120 cm	310-960-120
1,0 mm	160 cm	310-960-160
1,0 mm	180 cm	310-960-180
1,0 mm	230 cm	310-960-230
1,4 mm	120 cm	310-961-120
1,4 mm	160 cm	310-961-160
1,4 mm	180 cm	310-961-180
1,4 mm	230 cm	310-961-230
1,6 mm	120 cm	310-962-120
1,6 mm	160 cm	310-962-160
1,6 mm	180 cm	310-962-180
1,6 mm	230 cm	310-962-230
1,8 mm	120 cm	310-963-120
1,8 mm	160 cm	310-963-160
1,8 mm	180 cm	310-963-180
1,8 mm	230 cm	310-963-230
2,3 mm	120 cm	310-964-120
2,3 mm	160 cm	310-964-160
2,3 mm	180 cm	310-964-180
2,3 mm	230 cm	310-964-230
3,0mm	180cm	310-965-180
3,0mm	230cm	310-965-230

Standard delivery:

Metal tube distal tapered, metal handle.



On desire:

- ⇒ Metal tube Grounded, e.g. 310-000-230 S
- ⇒ Metal tube Isolated, e.g. 310-000-230 SI **NOTICE:** not for diameter 1.0 – 1.6 mm !
- ⇒ Plastic Handle, e.g. 310-000-230 P

ALL INSTRUMENTS ARE AUTOCLAVABLE
SPECIAL LENGTH ON REQUEST



Foreign Body Retriever, 2:2 teeth

Diameter	Length	Order-No.
1,8 mm	120 cm	310-970-120
1,8 mm	160 cm	310-970-160
1,8 mm	180 cm	310-970-180
1,8 mm	230 cm	310-970-230
2,3 mm	120 cm	310-971-120
2,3 mm	160 cm	310-971-160
2,3 mm	180 cm	310-971-180
2,3 mm	230 cm	310-971-230



Foreign Body Retriever, 2:1 teeth

Diameter	Length	Order-No.
1,8 mm	120 cm	310-972-120
1,8 mm	160 cm	310-972-160
1,8 mm	180 cm	310-972-180
1,8 mm	230 cm	310-972-230
2,3 mm	120 cm	310-973-120
2,3 mm	160 cm	310-973-160
2,3 mm	180 cm	310-973-180
2,3 mm	230 cm	310-973-230

Standard delivery:

Metal tube distal tapered, metal handle.



On desire:

- ⇒ Metal tube Grounded, e.g. 310-000-230 S
- ⇒ Metal tube Isolated, e.g. 310-000-230 SI **NOTICE:** not for diameter 1.0 – 1.6 mm !
- ⇒ Plastic Handle, e.g. 310-000-230 P

ALL INSTRUMENTS ARE AUTOCLAVABLE
SPECIAL LENGTH ON REQUEST



Grasping Forceps, serrated

Diameter	Length	Order-No.
1,8 mm	120 cm	310-974-120
1,8 mm	160 cm	310-974-160
1,8 mm	180 cm	310-974-180
1,8 mm	230 cm	310-974-230
2,3 mm	120 cm	310-975-120
2,3 mm	160 cm	310-975-160
2,3 mm	180 cm	310-975-180
2,3 mm	230 cm	310-975-230



Catheter Applying Forceps, 2,5mm

Diameter	Length	Order-No.
2,5 mm	120 cm	310-976-120
2,5 mm	160 cm	310-976-160
2,5 mm	180 cm	310-976-180
2,5 mm	230 cm	310-976-230

Standard delivery:

Metal tube distal tapered, metal handle.



On desire:

- ⇒ Metal tube Grounded, e.g. 310-000-230 S
- ⇒ Metal tube Isolated, e.g. 310-000-230 SI **NOTICE:** not for diameter 1.0 – 1.6 mm !
- ⇒ Plastic Handle, e.g. 310-000-230 P

ALL INSTRUMENTS ARE AUTOCLAVABLE
SPECIAL LENGTH ON REQUEST



Polycatch Forceps

Diameter	Length	Order-No.
2,5 mm	120 cm	310-977-120
2,5 mm	160 cm	310-977-160
2,5 mm	180 cm	310-977-180
2,5 mm	230 cm	310-977-230

Standard delivery:

Metal tube distal tapered, metal handle.



On desire:

- ⇒ Metal tube Grounded, e.g. 310-000-230 S
- ⇒ Metal tube Isolated, e.g. 310-000-230 SI **NOTICE:** not for diameter 1.0 – 1.6 mm !
- ⇒ Plastic Handle, e.g. 310-000-230 P

ALL INSTRUMENTS ARE AUTOCLAVABLE
SPECIAL LENGTH ON REQUEST



Three Prongs Retriever, sharp

Diameter	Length	Order-No.
1,0 mm	120 cm	310-978-120
1,0 mm	160 cm	310-978-160
1,0 mm	180 cm	310-978-180
1,0 mm	230 cm	310-978-230
1,8 mm	120 cm	310-979-120
1,8 mm	160 cm	310-979-160
1,8 mm	180 cm	310-979-180
1,8 mm	230 cm	310-979-230
2,3 mm	120 cm	310-980-120
2,3 mm	160 cm	310-980-160
2,3 mm	180 cm	310-980-180
2,3 mm	230 cm	310-980-230

Standard delivery:

Metal tube, metal handle



On desire:

- ⇒ Metal tube Grounded, e.g. 310-000-230 S
- ⇒ Spiral Isolated, e.g. 310-000-230 SI **NOTICE:** not for diameter 1.0 – 1.6 mm !
- ⇒ Plastic Handle, e.g. 310-000-230 P

ALL INSTRUMENTS ARE AUTOCLAVABLE
SPECIAL LENGTH ON REQUEST



Three Prongs Retriever, round

Diameter	Length	Order-No.
1,8 mm	120 cm	310-981-120
1,8 mm	160 cm	310-981-160
1,8 mm	180 cm	310-981-180
1,8 mm	230 cm	310-981-230
2,3 mm	120 cm	310-982-120
2,3 mm	160 cm	310-982-160
2,3 mm	180 cm	310-982-180
2,3 mm	230 cm	310-982-230



Four Prongs Retriever, sharp

Diameter	Length	Order-No.
2,3 mm	120 cm	310-983-120
2,3 mm	160 cm	310-983-160
2,3 mm	180 cm	310-983-180
2,3 mm	230 cm	310-983-230

Standard delivery:

Metal tube, metal handle



On desire:

- ⇒ Metal tube Grounded, e.g. 310-000-230 S
- ⇒ Spiral Isolated, e.g. 310-000-230 SI **NOTICE:** not for diameter 1.0 – 1.6 mm !
- ⇒ Plastic Handle, e.g. 310-000-230 P

ALL INSTRUMENTS ARE AUTOCLAVABLE
SPECIAL LENGTH ON REQUEST



Suture Cutting, stamp type

Diameter	Length	Order-No.
2,3 mm	120 cm	310-984-120
2,3 mm	160 cm	310-984-160
2,3 mm	180 cm	310-984-180
2,3 mm	230 cm	310-984-230



Suture Cutting, scissors type, 90° opening

Diameter	Length	Order-No.
2,3 mm	120 cm	310-985-120
2,3 mm	160 cm	310-985-160
2,3 mm	180 cm	310-985-180
2,3 mm	230 cm	310-985-230

Standard delivery:

Metal tube distal tapered, metal handle



On desire:

- ⇒ Metal tube Grounded, e.g. 310-000-230 S
- ⇒ Spiral Isolated, e.g. 310-000-230 SI **NOTICE:** not for diameter 1.0 – 1.6 mm !
- ⇒ Plastic Handle, e.g. 310-000-230 P

ALL INSTRUMENTS ARE AUTOCLAVABLE
SPECIAL LENGTH ON REQUEST



Hook scissor, single action

Diameter	Length	Order-No.
2,5 mm	120 cm	310-986-120
2,5 mm	160 cm	310-986-160
2,5 mm	180 cm	310-986-180
2,5 mm	230 cm	310-986-230

Standard delivery:

Metal tube distal tapered, metal handle



On desire:

- ⇒ Metal tube Grounded, e.g. 310-000-230 S
- ⇒ Spiral Isolated, e.g. 310-000-230 SI **NOTICE:** not for diameter 1.0 – 1.6 mm !
- ⇒ Plastic Handle, e.g. 310-000-230 P

ALL INSTRUMENTS ARE AUTOCLAVABLE
SPECIAL LENGTH ON REQUEST



3- Prongs Retriever, in Tefloncatheter for Universal Handle

Diameter	Length	Order-No.
2,3 mm	180 cm	310-987-180
2,3 mm	230 cm	310-987-230



4- Prongs Retriever, in Tefloncatheter for Universal Handle

Diameter	Length	Order-No.
2,3 mm	180 cm	310-988-180
2,3 mm	230 cm	310-988-230

UNIVERSAL HANDLE FOR LL-CONNECTOR, HUB 70 MM



Standard delivery:

- ⇒ In Tefloncatheter with Luer-Look Connector for Universal Handle
- ⇒ Without Handle

ALL INSTRUMENTS ARE AUTOCLAVABLE
SPECIAL LENGTH ON REQUEST



5- Prongs Retriever, in Tefloncatheter for Universal Handle

Diameter	Length	Order-No.
2,3 mm	180 cm	310-989-180
2,3 mm	230 cm	310-989-230



5- Prongs Retriever, in Tefloncatheter for Universal Handle

Diameter	Length	Order -No.
2,3 mm	180 cm	310-990-180
2,3 mm	230 cm	310-990-230

UNIVERSAL HANDLE FOR LL-CONNECTOR, HUB 70 MM



Standard delivery:

- ⇒ In Tefloncatheter with Luer-Look Connector for Universal Handle
- ⇒ Without Handle

ALL INSTRUMENTS ARE AUTOCLAVABLE
SPECIAL LENGTH ON REQUEST



3- Prongs Retriever, in Tefloncatheter for Universal Handle

Diameter	Length	Order-No.
2,3 mm	180 cm	310-991-180
2,3 mm	230 cm	310-991-230



4- Prongs Retriever, in Tefloncatheter for Universal Handle

Diameter	Length	Order-No.
2,3 mm	180 cm	310-992-180
2,3 mm	230 cm	310-992-230


UNIVERSAL HANDLE FOR LL-CONNECTOR, HUB 70 MM



Standard delivery:

- ⇒ In Tefloncatheter with Luer-Look Connector for Universal Handle
- ⇒ Without Handle

ALL INSTRUMENTS ARE AUTOCLAVABLE
SPECIAL LENGTH ON REQUEST



Four Prongs Retriever, in Tefloncatheter for Universal Handle

Diameter	Length	Order-No.
2,3 mm	120 cm	310-993-120
2,3 mm	160 cm	310-993-160
2,3 mm	180 cm	310-993-180
2,3 mm	230 cm	310-993-230

Standard delivery:



Metal tube distal tapered, metal handle



On desire:

- ⇒ Metal tube Grounded, e.g. 310-000-230 S
- ⇒ Spiral Isolated, e.g. 310-000-230 SI **NOTICE:** not for diameter 1.0 – 1.6 mm !
- ⇒ Plastic Handle, e.g. 310-000-230 P

ALL INSTRUMENTS ARE AUTOCLAVABLE
SPECIAL LENGTH ON REQUEST



Cleaning Brushes

311-500-102 for cleaning of suction and working channels of gastroscopes, 30mm, \varnothing 5mm, \varnothing metal-sheath2mm, 180cm

311-500-103 for cleaning of suction and working channels of colonoscopes, 30mm, \varnothing 5mm, \varnothing metal-sheath2mm, 230cm