

Unter Haßlen 20 · D-78532 Tuttlingen
Phone: +49 (0)7462/9470-0
Fax: +49 (0)7462/9470-50
e-mail: sales@rz-medizintechnik.com
<http://www.rz-medizintechnik.com>



Advanced RZ Endoscopic CLIP APPLYING



12/10

ADVANCED ENDOSCOPIC RZ CLIP APPLYING FORCEPS for laparoscopic and XXL Bariatric laparoscopic procedures



Advanced RZ Clip Applying Forceps, 360° rotatable, LL-rinsing

Article-No.	Working length	Diameter	clip size
300-803-001	330 mm / 13"	10 mm	● small
300-803-002	330 mm / 13"	10 mm	● medium
300-803-004	330 mm / 13"	10 mm	● medium - large
300-803-006	330 mm / 13"	12 mm	● large
300-803-401	450 mm / 17 ³ / ₄ "	10 mm	● small
300-803-402	450 mm / 17 ³ / ₄ "	10 mm	● medium
300-803-404	450 mm / 17 ³ / ₄ "	10 mm	● medium - large
300-803-406	450 mm / 17 ³ / ₄ "	12 mm	● large

- ADVANCED RZ Clip Applying Forceps take out RZ Clips easily.
- Clips are well positioned in the center of the jaws.
- Jaws hold clips securely in place, easy to clean through integrated flush channel.
- Light weight handle.





Curved tip design



Jaw for small clips



Jaw for medium and large clips



small
blue color code



medium
white color code



medium-large
green color code



large
yellow color code



300-400-901
Silicon cap for closing
luer-lock channel



- Clips does not move inside
- RZ clips closes securely and holds the vessel in the clip without rupture
- Grooves are non-traumatic

ENDOSCOPIC RZ CLIPS



| 300-880-100 Titanium Ligating Clip, Small size (S), package with 36 cartridges each 6 clips, blue colour code



| 300-880-200 Titanium Ligating Clip, Medium size (M), package with 36 cartridges each 6 clips, white colour code



| 300-880-300 Titanium Ligating Clip, Medium-Large size (ML), package with 18 cartridges each 6 clips, green colour code



| 300-880-400 Titanium Ligating Clip, Large size (L), package with 18 cartridges each 6 clips, yellow colour code

